

AFFINITY PURIFICATION OF SOLUBLE GST-TAGGED PROTEINS

www.ptglab.com

1. Lyse cells:

- Suspend the cell pellet in 30–35 ml of glutathione S-transferase (GST)-washing buffer with 10 mM PMSF and 0.5 M EDTA.
- b. Sonicate cells in an ice-bath at 200 W for 6 min.
- c. Rotate the lysed solution for 1 h at 4°C.
- d. Centrifuge the cell lysate for approximately 13 min at 8000 rpm, $4^{\circ}\text{C}.$

2. Bind protein to beads:

- a. Transfer the supernatant to $600 \mu l$ of GST-beads.
- b. Rotate the mixture overnight at 4°C.
- Collect the beads by centrifugation at 2000 rpm for 10–30 seconds, 4°C.
 Collect the protein bound beads in eppendorf tubes.

3. Wash out the unbound proteins from beads:

a. Wash the beads 3 times with 1ml of GST-washing buffer. Discard the supernatant.

4. Elute proteins from beads:

- a. Add 300 µl of GST-elution buffer to the beads.
- b. Rotate the mixture for 1 h at 4° C.
- c. Collect the supernatant by centrifugation at 300 rpm for 10-30 seconds.
- d. Repeat steps 4 a-c.
- e. Combine the eluent (total volume of 600 µl).
- f. Check the molecular weight and purity of the enriched protein by SDS-PAGE analysis.



AFFINITY PURIFICATION OF SOLUBLE GST-TAGGED PROTEINS

www.ptglab.com

Buffers Needed

GST-washing buffer (PBST buffer)	1000 ml
58 mM Na ₂ HPO ₄	8.24 g
17 mM NaH ₂ PO ₄	2.04 g
68 mM NaCl	3.98 g
1% Triton X-100	10 ml
Add ddH₂O to 1000 ml	
Adjust to pH 7.4	

GST-elution buffer	1000 ml
100 mM GSH	30.70 g
10% Glycerol	100 ml
1x PBST buffer	900 ml
Adjust to pH 8.0	