

WESTERN BLOTTING

Technical Tips and Troubleshooting

CONTENTS

| | |
|--------------|------------------------------------|
| <u>3</u> | Principles Of Western Blotting |
| <u>4-5</u> | Antigen Optimization |
| <u>6-8</u> | Low Molecular Weight Proteins |
| <u>9</u> | Blocking |
| <u>10-11</u> | Transfer |
| <u>12-13</u> | Membrane Choice |
| <u>14</u> | Monoclonal Or Polyclonal Antibody? |
| <u>15-16</u> | Controls |
| <u>17-20</u> | Troubleshooting |
| <u>21</u> | Contact Us |

PRINCIPLES OF WESTERN BLOTTING

- Western Blotting is an **Immunoblotting Technique**.
- Western Blotting is an analytical technique to detect **specific proteins** in a cell extract or tissue homogenate.
- Western Blotting relies on the specific binding between the protein-of-interest and an **antibody raised against this particular protein**.
- Main steps are:
 - Sample preparation**
 - Gel electrophoresis**
 - Transfer**
 - Blocking**
 - Detection**
 - Analysis**

ANTIGEN OPTIMIZATION

- Cell lysis is the breaking down of the cell membrane and the separation of proteins from the non-soluble parts of the cell.
- Dependent on the location of the protein of interest, a different lysis buffer is needed.

| Location | Lysis Buffer |
|--------------------------------------|-------------------------|
| Nuclear, Mitochondria | RIPA, Fraction protocol |
| Membrane bound, Whole cell lysate | RIPA or NP-40 |
| Cytoplasm | RIPA, Tris, HCL/Triton |

RIPA: 25mM Tris, HCl (pH 7.6) 150mM NaCl, 1% NP-40, 1% sodium deoxycholate, 0.1% SDS.


NP-40: 50 mM Tris, HCl (pH 8.5) 150 mM NaCl, 1% detergent.

ANTIGEN OPTIMIZATION

- The non-specific binding is favoured when a low amount of antigen is present
- For low abundance proteins the following applies:

Insufficient amount of the antigen might lead to non-specific binding. Load more protein.

Protein loading amount can help to optimize signal issues.




LOW MOLECULAR WEIGHT PROTEINS

Standard Gel:

- Glycin-Tris buffer system
- Detection of 30-250 kDa proteins

Gel for low molecular weight proteins:

- Tris-Tricine buffer
- Detection of <30 kDa proteins



H Schägger and G von Jagow. Tricine-sodium dodecyl sulfate-polyacrylamide gel electrophoresis for the separation of proteins in the range from 1 to 100 kDa. Anal Biochem. 1987;166(2):368-79.


LOW MOLECULAR WEIGHT PROTEINS

Tricine Gel Recipe For Low MW Proteins

| Items | Separation | Intermediate | Stacking |
|----------------------------|------------|--------------|----------|
| Gel Percentage | 15% | 10% | 4% |
| Gel Volume | 6 ml | 3 ml | 2 ml |
| 38% Glycerol | 1.6 | — | — |
| ddH ₂ O | — | 1.2 | 1.4 |
| 30% Acrylamide | 2.7 | 0.8 | 0.3 |
| 3.0 M Tris HCl (pH 8.5) | 2.14 | 1 | — |
| 1.0 M Tris HCl (pH 6.8) | — | — | 0.3 |
| 10% SDS | 0.06 | 0.03 | 0.02 |
| 10% APS | 0.06 | 0.03 | 0.02 |
| TEMED | 0.003 | 0.003 | 0.002 |

LOW MOLECULAR WEIGHT PROTEINS

Tricine Gel Recipe For Low MW Proteins




BLOCKING

- As the membrane has the ability to bind proteins, there might be some non-specific bindings.
- Blocking is necessary to optimize the signal-to noise ratio.
- **Bovine Serum Albumin (BSA)** or **non fat milk** with a low percentage of detergent are used for blocking.
- For phospho-proteins BSA is suggested
- Controls to optimize the best blocking reagent:
 - blocking buffer + substrate**
 - blocking buffer + 2.AB + substrate**
 - blocking buffer + 1.AB + 2.AB + substrate**

TRANSFER

- Proteins are transferred from the gel to a membrane to make them accessible for antibody detection.
- Factors that influence the quality of the protein transfer:
 - Transfer length**
 - Temperature**
 - Membrane**
 - Buffer composition**
- Membrane choice and pore size:
 - Pore size:**
Dependent on MW of target protein.
 - Membrane type:**
Dependent on the hydrophilic/ hydrophobic character of the target protein (*see membrane choice*).
- Low molecular weight proteins tend to “over transfer”:
Transfer slowly and reduce voltage and time.

TRANSFER



A representation of the components of a transfer "sandwich." Note the orientation of the gel and membrane: the PVDF membrane is situated nearest to the positive electrode. The binding of SDS to proteins results in the complex having an overall negative charge. Therefore, the SDS-bound proteins travel towards the positive electrode.




MEMBRANE CHOICE

- Protein binding is based on **hydrophobic interactions and charged interactions.**
- The **PVDF** and the **Nitrocellulose** membrane are the most common used membranes.

| PVDF (Polyviniliden Difluoride) | Nitrocellulose |
|---|---|
| Less fragile. | Fragile. |
| Soak in ethanol, isopropanol or methanol before use. | Soak in transfer buffer before use. |
| More expensive. | Less expensive. |
| Can be stripped and reused a lot of times. | Can't be reused as much. |
| Preferable for proteins with a hydrophilic character. | Preferable for proteins with a hydrophobic character. |

MEMBRANE CHOICE

PVDF






Beta-Catenin ●

DLMDGLPPGDSNQLAWFDTL

Works on both PVDF & NC

Proteins with a hydrophilic character work better on PVDF.

Nitrocellulose



NALP-3 ●

CLLQNLGLSEMYFNKETSA

Works better on NC

Proteins with a hydrophobic character work better on PVDF.

MONOCLONAL OR POLYCLONAL ANTIBODY

Monoclonal

- High specificity.
- Less background.
- Better reproducibility.
- Sensitive to epitope loss/changes.
- Batch consistency.



Polyclonal

- Higher affinity and comparable specificity.
- Better results with non-denatured samples.
- More tolerant of minor changes.
- Species transferable.
- Lot to lot variation.



CONTROLS

- Control samples on a Western Blot might help to understand the cause of non-specific bands.

Positive Control

- Cells overexpressing the target protein.
- Purified recombinant protein.
- Samples with proven positive signal.

Negative Control

- Samples targeted with RNA interference.
- Knockout samples.
- Samples with proven negative signal.

LOADING CONTROL ANTIBODIES

- Loading controls are used for semi-quantification of protein levels.
- They ensure that the observed alteration in target protein is due to experimental manipulations.
- Characteristics of a loading control: high expression, constitutive, unchanged expression during experiment.

Selected loading controls for different cell location

| Whole cell/ cytoplasmic | Nuclear | Mitochondrial |
|----------------------------|-----------------|------------------------|
| Actin, 42 kDa | Lamin B1,66 kDa | COX-4, 17 kDa |
| GAPDH, 36 kDa | PCNA,36 kDa | VDCA1/Porin, 31 kDa |
| Tubulin, 50-55 kDa | TBP, 38 kDa | – |

TROUBLESHOOTING

General Difficulties

| Potential cause | Suggested solution |
|--|--|
| Dark spots. | Antibodies binding to the blocking reagent/filter the blocking buffer. |
| Sample overloading. | Decrease the total protein loading. |
| Too much antibody. | Titrate the antibody concentration. |
| Inverse staining (i.e., white bands on a dark blot). | Too much primary and/or too much secondary antibody. |
| Uneven bands. | Uneven gel composition/insufficient buffer being added to the tank. |
| Blank areas/white spots. | Improper/uneven transfer. |
| “Smiling” bands. | Migration through the gel was too hot or too fast. |

TROUBLESHOOTING

Non-specific Binding

| Potential cause | Suggested solution |
|---|--|
| Sample degradation. | Prepare fresh lysates. Include protease inhibitors. |
| Interference from other isoforms. | Check literature for known isoforms. Use isoform-specific antibody. |
| Target protein abundance is lower than threshold of non-specific binding. | Load more protein. Enrich low abundance proteins by IP. |
| Inefficient SDS-PAGE separation. | Change the gel percentage to suit the target protein's MW. |

TROUBLESHOOTING

Weal/No Signal

| Potential cause | Suggested solution |
|--|--|
| Issues with the primary and/or secondary antibody. | Titrate the antibody. Change incubation time and temperature. Antibody may have lost activity. |
| Membrane choice. | Select PVDF or NC membranes based on hydrophobicity. Hydrophilicity of the target antigen. |
| Sodium azide contamination. | The presence of sodium azide inhibits the activity of HRP. Ensure sufficient washing. |
| Detection reagent not sensitive enough. | Dilute chemiluminescent reagents in high-purity water. Check several exposure times to achieve optimum detection. |

TROUBLESHOOTING

High Background

| Potential cause | Suggested solution |
|--|--|
| Inadequate washing. | Increase washing time and volume. |
| Dry membrane. | Ensure membrane does not dry out. |
| Film exposure too long/ detection reagent too sensitive. | Check different types and dilution of the detection reagent. |
| Insufficient blocking. | Increase the concentration of blocking reagent. Increase blocking time. Add Tween 20 to the blocking buffer. |
| Antibody concentration too high. | Use higher antibody dilution. |

CONTACT US

Proteintech Group
US Head Office

proteintech@ptglab.com

Proteintech Europe
United Kingdom

europe@ptglab.com

Proteintech
China Office

service@ptglab.com

Support

Available 24 hours via Live Chat and 9–5 (CDT) via phone.

Please visit us at www.ptglab.com for more information about our antibodies and technical tips.