

INDIRECT ELISA

www.ptglab.com

All steps are carried out at room temperature unless stated otherwise.
Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

-
- 1. Antigen coating:**
 - Dilute purified antigens to a final concentration of 0.2 µg/ml in **antigen-coating buffer** and add 100 µl of diluted antigen to each well of a 96-well ELISA plate.
 - Carefully cover the plate with adhesive plastic and incubate at 4°C overnight.

 - 2. Blocking:**
 - Empty the wells of **antigen-coating buffer** and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.
 - Add 200 µl **blocking buffer** per well to block residual protein-binding sites. Cover the plate with adhesive plastic and incubate for 1–2 h at 37°C.

 - 3. Antibody incubation:**
 - Dilute your primary antibody of choice with **blocking buffer** in a series e.g. 1:500, 1:1000, 1:2000, 1:4000 and so on, empty the wells of **blocking buffer** and then add 100 µl of each dilution per well. Repeat in duplicate, or triplicate, for accuracy. Cover the plate with adhesive plastic and incubate for 1 h at 37°C.
 - Empty the wells and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.
 - Dilute the HRP-conjugated secondary antibody with **blocking buffer** at an optimal concentration (a dilution factor within 1:10,000-1:100,000 is recommended) and add 100 µl of secondary antibody solution to each well. Cover the plate with adhesive plastic and incubate for 1 h at 37°C.
 - Empty the wells and wash 3 times with 200 µl **PBST buffer** for 5 minutes each time.

 - 4. Signal detection:**
 - Add 100 µl TMB substrate (mix equal volumes of **TMB buffer A** and **buffer B**) to each well with a multichannel pipette. Color development should peak after 15 minutes, at which time it should be stopped by adding 100 µl of 2 M H₂SO₄ per well. Read absorbance at 450 nm.

INDIRECT ELISA

www.ptglab.com

Buffers Needed

| Antigen coating-buffer | For 1000 ml | PBST buffer | For 1000 ml |
|-----------------------------------|-------------|---|-------------|
| 100 mM NaHCO ₃ | 8.4 g | 10 mM Na ₂ HPO ₄ | 1.42 g |
| Adjust pH to 9.6 | | 1.8 mM NaH ₂ PO ₄ | 0.22 g |
| Add ddH ₂ O to 1000 ml | | 140 mM NaCl | 8.19 g |
| | | 0.2 % Tween 20 | 2 ml |
| | | Adjust pH to 7.4 | |
| | | Add ddH ₂ O to 1000 ml | |

| Blocking buffer | For 100 ml |
|---------------------------|------------|
| 5% non-fat dry milk | 5 g |
| Add PBST buffer to 100 ml | |

| TMB buffer A | For 500 ml | TMB buffer B | For 500 ml |
|-----------------------------------|------------|------------------------------------|------------|
| NaAc•3H ₂ O | 13.6 g | TMB (first dissolved in 3 ml DMSO) | 0.15 g |
| Citric acid | 1.6 g | EDTA-2Na | 0.2 g |
| 30% H ₂ O ₂ | 0.3 ml | Citric acid | 0.95 g |
| Add ddH ₂ O to 500 ml | | Glycerol | 50 ml |
| | | Add ddH ₂ O to 500 ml | |