

## Flow Cytometry Triton X-100 Permeabilization Protocol

## Reagents required:

4% PFA Fix Buffer (4% paraformaldehyde dissolved in 1x PBS, adjust pH to 7.4)

0.1% Triton X-100 (in 1x PBS)

1x PBS

Flow cytometry antibodies

## **Experiment procedures:**

- 1. Harvest cells and wash them twice with 1x PBS by centrifugation at 350-500 x g for 5 minutes each time, discard the supernatant.
- 2. (Optional) Perform cell surface staining with recommended amount of fluorochrome-conjugated primary antibody, wash the cells with 1 mL staining buffer by centrifugation at  $350-500 \times g$  for 5 minutes, discard the supernatant.
- 3. Resuspend the cells with 200  $\mu$ L of 4% PFA Fix Buffer and vortex briefly, incubate for 20 minutes at room temperature in the dark.
- 4. Centrifuge at  $350-500 \times g$  for 5 minutes, discard the supernatant. Wash the cells twice with  $1 \times PBS$  by centrifugation at  $350-500 \times g$  for 5 minutes each time.
- 5. Resuspend the cells at a cell density of approximately 1 x 10<sup>6</sup> cells in 1 mL of 0.1% Triton X-100, incubate for 15 minutes at room temperature in the dark.
- 6. Centrifuge at 350-500 x g for 5 minutes, discard the supernatant. Wash the cells twice with 1x PBS by centrifugation at 350-500 x g for 5 minutes each time.
- 7. Resuspend the cells with 1x PBS.
- 8. Incubate the cells with the primary antibody in each 100  $\mu$ L of cell resuspension. The concentration of the primary antibody is based on the recommendations or the results of titration.
- 9. Incubate for 45-60 minutes at 4°C in the dark.
- 10. Wash the cells with 1 mL staining buffer by centrifugation at 350-500 x g for 5 minutes, discard the supernatant. Repeat.

Note: If using fluorochrome-conjugated primary antibodies, skip to step 14.

- 11. Resuspend the cells with diluted fluorochrome-conjugated secondary antibody in 100  $\mu$ L 1x PBS (use recommended concentration for secondary antibody dilution).
- 12. Incubate for 45-60 minutes at 4°C in the dark.
- 13. Wash the cells with 1x PBS by centrifugation at 350-500 x g for 5 minutes, discard the supernatant.
- 14. Resuspend the cells with 200-500  $\mu$ L 1x PBS and analyze on flow cytometer.