


GFP-Trap[®]

Immunoprecipitation of GFP-fusion proteins and their interacting factors with anti-GFP Nanobody conjugated to beads.

- Fast, reliable & efficient one-step immunoprecipitation
- For IP, Co-IP, MS & ChIP/RIP
- Ready-to-use
- No heavy & light antibody chains
- Compatible with almost all buffers
- More than 2,200 publications



Immunoprecipitation of GFP by GFP-Trap compared with a conventional anti-GFP antibody coupled to Protein A/G beads analyzed by SDS-Page. I: Input, FT: Flow-Through, B: Bound



**For more information & free sample
visit our website or contact:**

ChromoTek

E: order@chromotek.com

Proteintech Group

E: proteintech@ptglab.com



SCAN ME